

HAMA-BACTA™

مطهر و منظف الأرضيات
و الحوائط

Floor Cleaner & Broadsurface, Disinfectant Concentrate



Lemon Aroma

Concentrated disinfectant specially formulated for the cleaning and disinfection of floors and other large surfaces in hospitals

Features

Water-based disinfectant cleaner with Quaternary Ammonium Compounds as active ingredients. Contains surface active agents for its cleaning properties. Concentrated, aldehyde-free and phenol-free formulation Effective against bacteria, fungi and viruses (including HIV/HBV). Easily biodegradable Safe to use Action time of 5-60 minutes.

Directions for Use

10% solution: Add 100 ml of concentrate disinfectant to 1 litres of warm water and allow to mix. Using a floor type sponge cleaner, wipe all surfaces and allow to dry. No need to rinse surfaces can be allowed to dry naturally.

Composition

Aqua Water, Didecyldimethylammonium Chloride, Ethylenediaminetetraacetic Acid Nonionic Surfactant, Lemon Perfume and Other Substances.

Physical Properties

Green liquid, odour Lemon, pH 10.0 - 11.0. Soluble in water.

Packaging

1111-025 5 Litre Container

Standard Dilution

20 ml With 8 Litre from water
40 ml With 16 Litre from water
60 ml With 24 Litre from water
80 ml With 32 Litre from water
100 ml With 40 Litre from water

